



Product Datasheet

CD123 antibody (Biotin) (orb179857)

Description	Mouse monoclonal antibody conjugated to Biotin which recognizes CD123.
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, IHC-Fr, IP, WB
Immunogen	IL3 receptor alpha chain expressed on the surface of transiently transfected COS cells
Target	CD123
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-12 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P26951
Entrez	3563
Dilution Range	Flow cytometry: Recommended dilution: 1-12 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...



Flow
cytometry
analysis
(surface
stainin...