



Product Datasheet

CD161 antibody (orb179861)

Description	Mouse monoclonal antibody to CD161.
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, WB
Immunogen	human NK cells
Target	CD161
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q12918
Entrez	3820
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of
human
CD161
positive
lymph...



Flow
cytometry
surface
staining
pattern ...