



Product Datasheet

CD1c antibody (orb179867)


Description	Mouse monoclonal antibody to CD1c.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	human thymocytes
Target	CD1c
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P29017
Entrez	911
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human CD1c
positive
CD19 p...



Flow
cytometry
surface
staining
pattern ...