


Product Datasheet

CD271 antibody (Biotin) (orb179871)

Description	Mouse monoclonal antibody conjugated to Biotin which recognizes CD271.
Reactivity	Feline, Ferret, Human, Primate, Rabbit
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-P, IP, WB
Immunogen	Purified CD271 protein isolated from human melanoma cell line A875
Target	CD271
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Extracellular and intracellular staining; recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P08138
Entrez	4804
Dilution Range	Flow cytometry: Extracellular and intracellular staining; recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of REH cells from human periparturient...