


Product Datasheet

CD3 antibody (PerCP) (orb179874)

Description	Mouse monoclonal antibody against CD3 conjugated to PerCP
Reactivity	Mouse
Conjugation	PerCP
Tested Applications	FC
Immunogen	Mouse BM10-37 cytotoxic T lymphocytes
Target	CD3
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Isotype	Armenian Hamster IgG kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated Peridinin-Chlorophyll Protein (PerCP) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P22646
Entrez	12501
Dilution Range	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface stainin...