



Product Datasheet

CD45R antibody (FITC) (orb179889)

Description	Rat monoclonal antibody conjugated to FITC which recognizes CD45R.
Reactivity	Feline, Human, Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Abelson murine leukemia virus-induced pre-B tumor cells
Target	CD45R
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Rat IgG2a kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of murine splenocytes stained...



Flow cytometry surface staining pattern ...