


Product Datasheet

CD45R antibody (PE) (orb179891)

Description	Rat monoclonal antibody conjugated to PE which recognizes CD45R.
Reactivity	Feline, Human, Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	Abelson murine leukemia virus-induced pre-B tumor cells
Target	CD45R
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Rat IgG2a kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...