



Product Datasheet

CD79a antibody (FITC) (orb179902)


Description	Mouse monoclonal antibody conjugated to FITC which
Reactivity	Bovine, Equine, Gallus, Guinea pig, Human, Mouse, Opossum, Porcine, Rabbit, Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Synthetic peptide corresponding to amino acids 202-216 of human CD79a
Target	CD79a
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests. Intracellular staining.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P11912
Entrez	973
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests. Intracellular staining.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD79a
positive B
cel...



Flow
cytometry
intracellular
staining pa...