

Product Datasheet

CD79b antibody (FITC) (orb179903)

Description

Mouse monoclonal antibody conjugated to FITC which

Reactivity

Human

Conjugation

FITC

Tested

FC

Applications
Immunogen

Fraction of Ig-associated molecules isolated from Ramos B cells

Target

CD79b

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P40259

Entrez


974

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Expiration Date


12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD79b
positive B
cel...



Flow
cytometry
surface
staining
pattern ...