

Product Datasheet

CD79b antibody (PE) (orb179904)

Description

Mouse monoclonal antibody conjugated to PE which

Reactivity

Human

Conjugation

PE

Tested

FC

Applications
Immunogen

Fraction of Ig-associated molecules isolated from Ramos B cells

Target

CD79b

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P40259

Entrez


974

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date


12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human CD73
positive B
cell...



Flow
cytometry
surface
staining
pattern ...