



Product Datasheet

CD85j antibody (orb179908)


Description	Mouse monoclonal antibody to CD85j.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IP, WB
Immunogen	Hairy cell leukaemia cells
Target	CD85j
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q8NHL6
Entrez	10859
Dilution Range	Functional application: Blocking of HLA-G induced TGF-beta1 production. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow cytometry multicolor surface staini...



Separation of human CD85j positive B cel...



Flow cytometry surface staining pattern ...