



Product Datasheet

CLIC5a antibody (orb179913)

Description	Mouse monoclonal antibody to CLIC5a.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH-conjugated peptide corresponding to the amino acids 160-173 of human CLIC5a (Cys-Ile-Asp-Ala-Asn-Thr-Ser-Asp-Lys-Gly-Ser-Arg-Arg coupled with KLH).
Target	CLIC5
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9NZA1
Entrez	53405
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Western blotting analysis of CLIC5a in H...



Separation of NALM-6 cells stained using...