



Product Datasheet

HLA-A antibody (orb179915)

Description	Mouse monoclonal antibody to HLA-A
Reactivity	Bovine, Feline, Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC-Fr, IP, MC
Immunogen	Membrane of human tonsil cells
Target	HLA Class I
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Non-reducing conditions. Flow cytometry: Recommended dilution: 1-12 µg/ml
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-12 µg/ml
Expiration Date	12 months from date of receipt.



Separation of leukocytes stained using a...



Flow cytometry surface staining pattern ...