



Product Datasheet

IFN gamma antibody (Biotin) (orb179917)


Description	Mouse monoclonal antibody against IFN gamma conjugated to Biotin
Reactivity	Human, Primate
Conjugation	Biotin
Tested Applications	ELISA, FC, ICC, IHC-P, IP, WB
Immunogen	Interferon gamma derived from human leukocytes
Target	IFN gamma
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P01579
Entrez	3458
Dilution Range	Flow cytometry: Recommended dilution 1-4 µg/ml. Intracellular staining. ELISA: Detection antibody in combination with capture antibody NIB42.
Expiration Date	12 months from date of receipt.



Cytokine IFN gamma (IFN γ) sandwich ELIS...



Flow cytometry multicolor intracellular ...



Separation of human IFN-gamma positive C...