

Product Datasheet

Ig Lambda Light Chain Antibody (APC) (orb179920)

Description

Anti-Hu Ig Lambda Light Chain APC

Reactivity

Human

Conjugation

APC

Tested Applications

FC

Immunogen

Human IgA1 lambda myeloma protein

Target

Ig Lambda Light Chain

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Extracellular and intracellular staining.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity


Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Extracellular and intracellular staining.

Expiration Date

12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...