


Product Datasheet

NCK1 antibody (orb179922)

Description	Mouse monoclonal antibody to NCK1.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Purified recombinant second SH3 domain of human Nck1 fused to GST
Target	NCK1
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution: 2 µg/ml; positive control: Jurkat cell lysate. Immunocytochemistry: Recommended dilution: 2 µg/ml; positive control: human NCK1-transfected COS-7 cells.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P16333
Entrez	4690
Dilution Range	Western blotting: Recommended dilution: 2 µg/ml; positive control: Jurkat cell lysate. Immunocytochemistry: Recommended dilution: 2 µg/ml; positive control: human NCK1-transfected COS-7 cells.
Expiration Date	12 months from date of receipt.



Western blotting analysis of NCK1 in Jurkat cell lysate.