

Product Datasheet

Thyrotropin antibody (orb179927)

Description	Mouse monoclonal antibody to Thyrotropin.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, RIA
Immunogen	Human thyrotropin.
Target	hTSH
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	RIA: The antibody TSH-51 is suitable in combination with the antibody TSH-116 for immunometric assays in the screening of neonatal hypothyroidism.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	RIA: The antibody TSH-51 is suitable in combination with the antibody TSH-116 for immunometric assays in the screening of neonatal hypothyroidism.
Expiration Date	12 months from date of receipt.