

Product Datasheet

TNF alpha antibody (APC) (orb179928)

Description

Mouse monoclonal antibody against TNF alpha



Reactivity

Human, Porcine, Primate

Conjugation

APC

Tested

FC

Applications

Immunogen

Recombinant human TNF-alpha

Target

TNF alpha

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[P01375](#)

Entrez


[7124](#)

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Expiration Date

12 months from date of receipt.



Flow
cytometry
analysis
(intracellular
S...