

## Product Datasheet

### TNF alpha antibody (FITC) (orb179929)

**Description**

Mouse monoclonal antibody against TNF alpha

**Reactivity**

Human, Porcine, Primate

**Conjugation**

FITC

**Tested**

FC

**Applications**

**Immunogen**

Recombinant human TNF-alpha

**Target**

TNF alpha

**Preservatives**

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

**Storage**

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

**Note**

For research use only

**Application notes**

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests. Intracellular staining.

**Isotype**

Mouse IgG1 kappa

**Clonality**

Monoclonal

**Purity**

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

**Uniprot ID**

**P01375**

**Entrez**


**7124**

**Dilution Range**

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests. Intracellular staining.


**Expiration Date**

12 months from date of receipt.




Flow  
cytometry  
multicolor  
intracellular

...



Separation  
of human  
TNF alpha  
positive C...



Flow  
cytometry  
intracellular  
staining pa...