

Product Datasheet

Human ADAMDEC1 protein (orb180165)

Description	Recombinant human ADAMDEC1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,1mM ZnCl2,pH7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Lys206-Glu470
Source	Human Cells
Uniprot ID	O15204
Expiration Date	6 months from date of receipt.