

Product Datasheet

Human NECAP2 protein (orb180166)

Description	Recombinant human NECAP2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB ₁ 150mM NaCl,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS,2mM EDTA,1mM DTT,pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His&C-His
Note	For research use only
Application notes	6His tag at the N-terminus, 6His tag at the C-terminus.
Protein Sequence	Met1-Phe263
Source	E.coli
Uniprot ID	Q9NVZ3
Expiration Date	6 months from date of receipt.