

Product Datasheet

Human APOC2 protein (orb180173)

Description	Recombinant human APOC2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mMPB,150mM NaCL,30%glycerol,pH7.2.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS,50%glycerol,PH7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Glu82
Source	E.coli
Uniprot ID	P35354
Expiration Date	6 months from date of receipt.