

Product Datasheet

Human ABCB5 protein (orb180175)

Description	Recombinant human ABCB5 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB□150mM NaCl,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-Trx
Note	For research use only
Application notes	Trx tag at the N-terminus.
Protein Sequence	Ile141-Val247
Source	E.coli
Uniprot ID	Q2M3G0
Expiration Date	6 months from date of receipt.