

Product Datasheet

Human CAMK1D protein (orb180197)

Description	Recombinant human CAMK1D protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB□150mM NaCl,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-GST
Note	For research use only
Application notes	GST tag at the N-terminus.
Protein Sequence	Met1-Lys385
Source	E.coli
Uniprot ID	Q8IU85
Expiration Date	6 months from date of receipt.