

Product Datasheet

Human CSNK1G2 protein (orb180206)

Description	Recombinant human CSNK1G2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,500mM NaCl,10%Glycerol,1mM DTT, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 500mM NaCl, 10% Glycerol, 1mM DTT, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met18-Lys415
Source	E.coli
Uniprot ID	P78368
Expiration Date	6 months from date of receipt.