

Product Datasheet

CXCR2 Antibody (orb18026)

Description

CXCR2 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human CXCR2 (196-212aa CYEDMGNNNTANWRMLLR), different from the related mouse sequence by six amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

40759 MW

Uniprot ID


P25025

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Rat

Expiration Date

12 months from date of receipt.



Western Blot analysis of human carcinoma...