

Product Datasheet

Human Glucagon protein (orb180279)

Description	Human Glucagon protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, 50% Glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, 50% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Arg21-Lys180
Source	Human Cells
Uniprot ID	P01275
Expiration Date	6 months from date of receipt.