

Product Datasheet

Human HSPB7 protein (orb180284)

Description	Recombinant human HSPB7 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,0.2M NaCl,2mM DTT,50%glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 0.2M NaCl, 2mM DTT, 50%Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Ile170
Source	E.coli
Uniprot ID	Q9UBY9
Expiration Date	6 months from date of receipt.