

Product Datasheet

Human KHK protein (orb180325)

Description	Recombinant human KHK protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 50nM KCl, 10%Glycerol, pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 50mM KCl, 10% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Val298
Source	Human Cells
Uniprot ID	AAH06233.1
Expiration Date	6 months from date of receipt.