

## Product Datasheet

### Human LRG1 protein (orb180328)

|                          |                                                                                            |
|--------------------------|--------------------------------------------------------------------------------------------|
| <b>Description</b>       | Recombinant human LRG1 protein                                                             |
| <b>Endotoxins</b>        | Less than 0.1 ng/ug (1 EU/ug)                                                              |
| <b>Preservatives</b>     | Supplied as a 0.2 µm filtered solution of 20mM Tris,200mM NaCl,10%Glycerol,pH8.0.          |
| <b>Form/Appearance</b>   | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 200mM NaCl, 10% Glycerol, pH 8.0. |
| <b>Storage</b>           | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.     |
| <b>Tag</b>               | C-Fc&His                                                                                   |
| <b>Note</b>              | For research use only                                                                      |
| <b>Application notes</b> | Fc, 6His tag at the C-terminus.                                                            |
| <b>Protein Sequence</b>  | Val36-Gln347                                                                               |
| <b>Source</b>            | Human Cells                                                                                |
| <b>Uniprot ID</b>        | AAH70198.1                                                                                 |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                             |