

Product Datasheet

Human MIF protein (orb180337)

Description	Recombinant human MIF protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB,150mMNaCl,PH7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Ala115
Source	E.coli
Uniprot ID	P14174
Expiration Date	6 months from date of receipt.