

## Product Datasheet

### Human UBE2M protein (orb180352)

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Recombinant human UBE2M protein                                                                  |
| <b>Endotoxins</b>       | Less than 0.1 ng/ug (1 EU/ug)                                                                    |
| <b>Preservatives</b>    | Supplied as a 0.2 µm filtered solution of 50mM HEPES,pH 7.5,2mM DTT,150mM NaCl,10%glycerol .     |
| <b>Form/Appearance</b>  | Supplied as a 0.2 um filtered solution of 50mM HEPES, pH 7.5, 2mM DTT, 150mM NaCl, 10% Glycerol. |
| <b>Storage</b>          | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.           |
| <b>Tag</b>              | tag free                                                                                         |
| <b>Note</b>             | For research use only                                                                            |
| <b>Protein Sequence</b> | Met1-Lys183                                                                                      |
| <b>Source</b>           | E.coli                                                                                           |
| <b>Uniprot ID</b>       | <b>P61081</b>                                                                                    |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                   |