

Product Datasheet

Human LCN2 protein (orb180358)

Description	Recombinant human LCN2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS, 50% glycerol, pH 7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, 50% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Gln21-Gly198
Source	E.coli
Uniprot ID	P80188
Expiration Date	6 months from date of receipt.