

Product Datasheet

Human FKBP7 protein (orb180370)

Description	Recombinant human FKBP7 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM GaCl ₂ , 10% Glycerol, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM CaCl ₂ , 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Gln24-Leu222
Source	Human Cells
Uniprot ID	Q9Y680
Expiration Date	6 months from date of receipt.