

## Product Datasheet

### Human PDCD1LG2 protein (orb180380)

|                          |                                                                                                                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Recombinant human PDCD1LG2 protein                                                                                                                                                                                                    |
| <b>Endotoxins</b>        | Less than 0.1 ng/ug (1 EU/ug)                                                                                                                                                                                                         |
| <b>Preservatives</b>     | Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.                                                                                                                                                              |
| <b>Form/Appearance</b>   | Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.                                                                                                                                                           |
| <b>Storage</b>           | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| <b>Tag</b>               | C-His                                                                                                                                                                                                                                 |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                 |
| <b>Application notes</b> | 6His tag at the C-terminus.                                                                                                                                                                                                           |
| <b>Protein Sequence</b>  | Leu20-Pro219                                                                                                                                                                                                                          |
| <b>Source</b>            | Human Cells                                                                                                                                                                                                                           |
| <b>Uniprot ID</b>        | <b>Q9BQ51</b>                                                                                                                                                                                                                         |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                                        |