

Product Datasheet

Human RCN3 protein (orb180398)

Description	Recombinant human RCN3 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, 10%Glycerol, pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM DTT, 10% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Lys21-Leu328
Source	Human Cells
Uniprot ID	Q96D15
Expiration Date	6 months from date of receipt.