

Product Datasheet

Human ECP protein (orb180400)

Description	Human ECP protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,1mMDTT,10%Glycerol,pH7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mMDTT, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Arg28-Ile160
Source	Human Cells
Uniprot ID	AAH96060.1
Expiration Date	6 months from date of receipt.