

Product Datasheet

Human ACAT2 protein (orb180415)

Description	Recombinant human ACAT2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB ₁ 150mM NaCl,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl,100mM NaCl,1mM DTT,10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-Trx&His
Note	For research use only
Application notes	Trx, 6His tag at the N-terminus.
Protein Sequence	Met1-Glu397
Source	E.coli
Uniprot ID	O75908
Expiration Date	6 months from date of receipt.