

Product Datasheet

Human UBE2D3 protein (orb180452)

Description	Recombinant human UBE2D3 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT,10% glycerol, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl,150mM NaCl,10% Glycerol,pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Met147
Source	E.coli
Uniprot ID	AAH66917
Expiration Date	6 months from date of receipt.