

Product Datasheet

Human UBA5 protein (orb180454)

Description	Recombinant human UBA5 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl pH 8.0, 1mM DTT, 10% glycerol, 50mM NaCl.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 50mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Met404
Source	E.coli
Uniprot ID	Q9GZZ9
Expiration Date	6 months from date of receipt.