

## Product Datasheet

### CD206 antibody (orb180464)

**Description**

CD206 antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Porcine, Rat

**Conjugation**

Unconjugated

**Tested**

ELISA, IHC-P

**Applications**
**Immunogen**

KLH conjugated synthetic peptide derived from human CD206. Please contact us for the exact immunogen sequence. The peptide is available as orb374819.

**Target**

CD206

**Form/Appearance**

10 mM PBS, 0.02% sodium azide

**Concentration**

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

**Storage**

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Polyclonal antibodies are purified by peptide affinity chromatography

**MW**

166 kDa

**Uniprot ID**
**Q61830, P22897**
**NCBI**
**002429, 31, 002438**
**Entrez**
**4360**
**Dilution Range**


IHC-P: 1:100-500

**Expiration Date**


12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P staining of pig large intestines t...



Immunohistochemical staining of mouse sk...