

Product Datasheet

SNAIL antibody (orb180466)

Description

SNAIL antibody



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Western blot analysis
of Human breast ca...

Immunogen

KLH conjugated synthetic peptide derived from human SNAIL. Please contact us for the exact immunogen sequence. The peptide is available as orb374820.

Target

SNAIL

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

29 kDa

Uniprot ID

[O95863](#)

NCBI

[005976](#), [005985](#), [31](#)

Entrez


[6615](#)

Dilution Range


WB: 1:200, IHC-P: 1:500

Expiration Date


12 months from date of receipt.



Western blot analysis
of Human breast ca...



IHC-P staining of
human thyroid
carcinom...



Immunohistochemical
staining of paraffin...