


Product Datasheet

ATXN3 antibody (orb180467)

Description	Goat polyclonal antibody to ATXN3. This protein appears to be a component of the ubiquitin proteasom...
Species/Host	Goat
Reactivity	Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant peptide derived from within residues 120 aa to 250 aa of human ATXN3 produced in E. coli.
Target	ATXN3
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 or 500 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:250-1:2,000
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y cell lysate.