



Product Datasheet

Clathrin HC antibody (orb180469)

Description	Goat polyclonal to Clathrin Heavy Chain. Clathrin is formed by three clathrin heavy chains and three...
Species/Host	Goat
Reactivity	Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Purified recombinant peptide derived from within residues 50 aa to N-terminus of human CLTC produced in E. coli.
Target	Clathrin HC
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 µl (2 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	2 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:250-1:1,000, IF:1:25-1:250
Expiration Date	12 months from date of receipt.



Western blot analysis of At-T20 cell lin...



Confocal immunofluorescence analysis of ...