


Product Datasheet

GroEL antibody (orb180470)

Description	Goat polyclonal antibody to GroEL. GroEL protein is a 60 kDa chaperone that promotes refolding of mi...
Species/Host	Goat
Reactivity	Bacteria, Pseudomonas sp
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant peptide derived from within residues 400 aa to the C-terminus of Burkholderia GroEL produced in E. coli.
Target	GroEL
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:5,000
Expiration Date	12 months from date of receipt.



Western blot analysis of MBP-GroEL C-ter...