

Product Datasheet

SARS-CoV-2 (Omicron BA.2.75) Spike S1, Fc Tag (orb1805227)

Description	SARS-CoV-2 (Omicron BA.2.75) Spike S1, Fc Tag is expressed from HEK293. It contains AA Ser 13-Arg 68...
Species/Host	Human
Reactivity	Virus
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Spike S1
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Ser 13 - Arg 682
Isotype	Fc
Purity	> 85%
MW	202.64 kDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.