



Product Datasheet

Anti-Rabies virus GP Reference Antibody (orb1805415)

Description	Anti-Rabies virus GP Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and t...
Reactivity	Virus
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	Rabies virus GP
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-Rabies
virus GP
Reference
Antibody ...



The purity
of Anti-
Rabies
virus GP
Refer...