



Product Datasheet

Anti-CAIX Reference Antibody (orb1805432)

Description	Anti-CAIX Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light ch...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	CAIX
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	144.62 kDa
Expiration Date	12 months from date of receipt.



Anti-CA9 /
CAIX
Reference
Antibody
(gire...



The purity
of Anti-CA9 /
CAIX
Reference
...