



Product Datasheet

Anti-GPNMB Reference Antibody (orb1805465)

Description	Anti-GPNMB Reference Antibodyd antibody and Monomethyl auristatin E (MMAE) used for the study of bre...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	GPNMB
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-GPNMB
Reference
Antibody
(glembatum...



The purity of
Anti-GPNMB
Reference
Antib...