

Product Datasheet

Anti-HGFR / c-Met Reference Antibody (orb1805476)

Description

Anti-HGFR / c-Met Reference Antibodydy and

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

HGFR / c-Met

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%


MW

145.64 kDa


Expiration Date

12 months from date of receipt.


kDa NR M R



Anti-HGFR /
c-Met
Reference
Antibody
(te...



The purity
of Anti-
HGFR / c-
Met
Referenc...



Immobilized
human
cMet His at
2 µg/mL c...